Introcan Safety® Deep Access

Cannulation under ultrasound





Preparation & pre-scan

- Confirm all items are accessible for an ultrasound guided peripheral IVC insertion and device stabilization
- Align display of ultrasound device in direction of cannulation for a optimized hand-eye coordination
- Apply gel
- Identify appropriate vessel and insertion site


Pre-insertion

- Use an aseptic non-touch technique and disinfect puncture site according to organizational policies
- Cover probe & apply gel
- Find predefined vessel




Insert and observe

- Adjust angle of insertion and access vessel
- Track your needle tip through the tissue and confirm vessel access via ultrasound



Lower and advance slightly

- Lower catheter to 15° angle 1
- Advance the catheter and needle unit slightly to ensure catheter tip is in vessel
- Advance catheter • Use the push-off plate to advance catheter off the needle into the vessel





Occlude and remove needle

- Make sure to put down the ultrasound probe and free your hand to occlude the vein
- Remove the needle while occluding the vein and discard the shielded needle

Catheter Stabilization

- Connect a primed extension set to the catheter hub
- Check patency via aspiration and flush the catheter
- Stabilize catheter and apply sterile dressing according to organizational policies

B. Braun Melsungen AG | Hospital Care | 34209 Melsungen | Germany Tel. +49 5661 71-0 | www.bbraun.com

Date of last revision: 03.2023