

B BRAUN Hospital Care	<p align="center">List LP0001</p> <p align="center">Discofix® C</p>	Dokument-Nr.: LP001 Änderungsstand: 2 Seite: 1 von 1
-----------------------------------	---	--

In the IV-therapy are the pharmaceutical drugs below common.
All this indexed drugs have a significant potential to generate cracks in plastic materials. Discofix C can be used up to 96 hours with this substances without cracks.

Pharmaceutical Drug	Geneneric Name	Indication	Supplier
Phenhydan	Phenytoin-Natrium α-Tetrahydrofurfuryl-ω-hydroxyoligo-(oxethylen)-1(2)(Glycofurol 75), Trometamol (THAM), Editin-acid Dinatriumsalt, Natriumhydroxid	Hydantoine Antiepileptikum Antiarrhythmikum	Desitin
Sandimmun-Infusion-Solution	Ciclosporin Ethanol, Poly(oxyethylen)-40-Rizinusöl 33 Vol.% Alcohol	Immunsuppressivum	Sandoz
Intralipid / Lipo-fundin etc.	Lipide	Parenteral Nutrition	Kabi, B.Braun, usw.
Nimotop	Nimodipin Ethanol 96%. Macrogol 400, Natriumcitrat 2H2O, waterfree lemon-acid, 23,7 Vol.% Alcohol	Ischaemia	Desitin
Epanutin/Epinat	Phenytoin-Natrium Ethanol, Natriumhydroxid, Propylenglycol, e 10,7 Vol.% Alcohol	status epilepticus	Parke-Davis
Kodan-Spray	2-propanol, 1-Propanol, 2-Biphenylol, Hydrogen-Peroxide, Water	Skin disinfection	Schülke+ Mayr
Disoprivan 1 and 2% Propofol Lipuro etc.	Propofol	Anesthetic	Zeneca B.Braun
Vepesid J	Etoposid Benzylalkohol Macrogol 300 Ethanol absolut	Cytostatic	Bristol-Myers Squibb
Endoxan	Cyclophosphamid	Cytostatic	Asta Medica

	Name	Abteilung/Funktion	Datum	Unterschrift
erstellt:	Heil	Marketing and Development	24.09.1996	
geprüft:	Zihlmann	Development	24.09.1996	
Revised	Heil	Strategic Marketing	23.05.2003	