



The Medical Usability Sign for outstanding usability of a medical product was awarded on the 18th of September 2008 to:

Product / Category **OrthoPilot FS101 / Product and Hardware Design**
Manufacturer **Aesculap AG**

USE-LAB reasons why intended users can rely on an optimized ease of use of **Aesculap's OrthoPilot**:

- ✓ Packaging ensures effective protection during shipment.
- ✓ Transport position allows compact dimensions.
- ✓ Effective locking mechanism of the opto-electronical device ensures safe transport.
- ✓ Smooth-running, big wheels and comparatively low weight support easy handling when manoeuvring the device.
- ✓ Effective brakes ensure safe and stable stand.
- ✓ Circular handles in comfortable height offer good alternatives to take hold of the device.
- ✓ Three joints offer a good adjustability of the camera bar, especially in horizontal plane.
- ✓ OrthoPilot is easy to move and to manoeuvre.
- ✓ Status LEDs show the device's state distinctly.
- ✓ Height adjustable monitor
- ✓ The foot switch panel and the power cable can be easily and safely stored at the device.
- ✓ Plane surfaces of the device's housing allow easy cleaning.
- ✓ The chassis offers compact dimensions and stability.



MEDICAL USABILITY SIGN

