

PRODUCT DATA SHEET

Flexima®/D VZ R Key Convex Wafers Two-piece stoma bag systems



Flexima® Key is a flexible adhesive two-piece system with a guided coupling system for precise, simple and secure bag positioning on the wafer. The soft material of the coupling system provides flexibility which allows the wafer to adapt well to different body contours. Flexima® Key convex wafers are designed for flush or slightly retracted stomas to be used with the Flexima® Key closed, drainable or urostomy bags.

Features

- Gentle convexity of 4 mm
- Adhesive flexible coupling ring with guiding system allowing three different bag positions
- Slots for belt fixation, three different positions available
- Available in:
 - Cut-to-fit:
 - Ø50 mm: cut-to-fit Ø15–25 mm
 - Ø60 mm: cut-to-fit Ø15–35 mm and Ø15–45 mm
 - Pre-cut:
 - Ø60 mm: pre-cut Ø25 mm, Ø30 mm, Ø35 mm, Ø40 mm
 - Boxes of 5 wafers

Article Description	Article Number
SOFTIMA KEY CONVEX B PLATE 50/15-25 MM	63052N
SOFTIMA KEY CONVEX B PLATE 60/15-35 MM	63062N
SOFTIMA KEY CONVEX B PLATE 60/15-45 MM	63063N
SOFTIMA KEY CONVEX B PLATE 60/25 MM	63064N
SOFTIMA KEY CONVEX B PLATE 60/30 MM	63066N
SOFTIMA KEY CONVEX B PLATE 60/35 MM	63067N
SOFTIMA KEY CONVEX B PLATE 60/40 MM	63069N