

# PRODUCT DATA SHEET

## Flexima® Key Colostomy Bags 2-piece stoma bag systems



Flexima® Key is a flexible adhesive two-piece system with a guided coupling system for precise, simple and secure bag positioning on the wafer. The soft material of the coupling system provides flexibility which allows the wafer to adapt well to different body contours. The colostomy bags are ergonomically shaped and optimized in size.

### Features

#### Colors and capacity (ISO 8670-2):

- Ø40 mm: capacity of +/- 440 ml, beige and transparent
- Ø50 mm: capacity of +/- 420 ml, beige and transparent
- Ø60 mm: capacity of +/- 550 ml, beige and transparent
- Ø80 mm: capacity of +/- 520 ml, beige and transparent

#### Ergonomically shaped colostomy bags:

- Integrated laminar technology (LT) filter: optimized degassing and deodorizing capacity
- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)
- To be used with the Flexima® Key wafer corresponding to the chosen diameter of the bag
- Boxes of 30 bags

**Easy bag positioning and removal:**

- Guided positioning of the bag on the wafer thanks to the "key" system: The tab of the adhesive coupling ring on the bag is inserted into the slot of the wafer to guarantee a precise positioning
- A special cut-away section at the top of the flange facilitates the bag removal
- When the bag is removed no residue of the adhesive bag ring is left on the wafer

**Usage**

Ostomy appliance

Article Description	Article Number
SOFTIMA KEY CLOSED POUCH BEIGE 40 MM	63140N
SOFTIMA KEY CLOSED POUCH BEIGE 50 MM	63150N
SOFTIMA KEY CLOSED POUCH BEIGE 60 MM	63160N
SOFTIMA KEY CLOSED POUCH BEIGE 80 MM	63180N