Safety Data Sheet

according to UK REACH Regulation

Issue date: 22/08/2022 Revision date: 22/08/2022 Supersedes: 29/06/2020 Version: 1.1

SDS No: 00056-0401



SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product name : Prontosan Wound Gel X

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Use of the substance/mixture : For cleansing, decolonization and moisturizing of acute/chronic skin wounds and burns

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Manufacturer Supplier

B. Braun Medical AG B. Braun Medical Ltd.

Seesatz 17 Thorncliffe Park Estate, Brookdale Road

CH-6204 Sempach - Switzerland Sheffield S35 2PW, UK - Medical Information Department

T +41 (0) 58 / 258 50 00 T +44 (0)114 225 9000 <u>info.bbmch@bbraun.com</u> medinfo.bbmuk@bbraun.com

E-mail address of competent person responsible for the SDS: sds@gbk-ingelheim.de

1.4. Emergency telephone number

Emergency number : INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a)

In England and Wales: NHS 111 In Scotland: NHS 24 - dial 111

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to GB CLP Regulation

Hazardous to the aquatic environment – Chronic Hazard, Category 3 H412

Full text of H- and EUH-statements: see section 16

Adverse physicochemical, human health and environmental effects

Harmful to aquatic life with long lasting effects.

2.2. Label elements

Labelling according to GB CLP Regulation

Signal word (CLP) : -

Hazard statements (CLP) : H412 - Harmful to aquatic life with long lasting effects.

Precautionary statements (CLP) : P273 - Avoid release to the environment.

P501 - Dispose of contents and container to an approved waste disposal plant.

EUH-statements : EUH208 - Contains polyhexamethylene biguanide hydrochloride; PHMB. May produce an

allergic reaction.

Extra phrases : This product is not subject to labeling obligation according to art. 1 (5) d) of Regulation (EC)

no. 1272/2008 under the form of a medical product because it is used under body contact.

Labelling according to: exemption for packages of a capacity of 125ml or less

Hazard statements (CLP) : H412 - Harmful to aquatic life with long lasting effects.

EUH-statements : EUH208 - Contains polyhexamethylene biguanide hydrochloride; PHMB. May produce an

allergic reaction.

Extra phrases : This product is not subject to labeling obligation according to art. 1 (5) d) of Regulation (EC)

no. 1272/2008 under the form of a medical product because it is used under body contact.

2.3. Other hazards

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401



This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII Contains no PBT/vPvB substances ≥ 0.1% assessed in accordance with REACH Annex XIII

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Comments

: Mixture of the substances listed below with nonhazardous additives

Name	Product identifier		Classification according to GB CLP Regulation
polyhexamethylene biguanide hydrochloride; PHMB	(CAS-No.) 27083-27-8 (EC Index-No.) 616-207-00-X	0,1	Carc. 2, H351 Acute Tox. 2 (Inhalation), H330 Acute Tox. 4 (Oral), H302 STOT RE 1, H372 Eye Dam. 1, H318 Skin Sens. 1B, H317 Aquatic Acute 1, H400 (M=10) Aquatic Chronic 1, H410 (M=10)

Full text of H- and EUH-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general

: Data of item 4 do partly not refer to the use and the regular employing of the product (in this sense consult package leaflet and expert information), but to liberation of major amounts in case of accidents and irregularities. Take off immediately all contaminated clothing. Wash contaminated clothing before reuse.

First-aid measures after inhalation

: Remove person to fresh air and keep comfortable for breathing. In the event of symptoms refer for medical treatment.

First-aid measures after skin contact

: Get medical advice if skin irritation persists.

First-aid measures after eye contact

: Rinse thoroughly and plentifully with water, also under the eyelids. If eye irritation persists, consult a specialist.

First-aid measures after ingestion

: Rinse mouth. Drink plenty of water. Do not give an unconscious person anything to drink. Call a physician immediately. Do not induce vomiting without medical advice.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact

: May cause sensitisation of susceptible persons by skin contact.

Symptoms/effects after eye contact

: May cause eye-irritation of susceptible persons.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Product does not burn, fire-extinguishing activities according to surrounding.

Unsuitable extinguishing media : Do not use a solid water stream as it may scatter and spread fire.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Non flammable.

Printing date: 28/06/2023 GB - en 2/10

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401

Explosion hazard : Product is not explosive.

Hazardous decomposition products in case of fire : Toxic fumes may be released.

5.3. Advice for firefighters

Precautionary measures fire : Cool endangered containers with water spray jet.

Firefighting instructions : Move containers from fire area if it can be done without personal risk. Use water spray or

fog for cooling exposed containers.

Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained

breathing apparatus. Complete protective clothing.

B BRAUN

Other information : Fire residues and contaminated firefighting water must be disposed of in accordance with

the local regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Ensure adequate ventilation. Evacuate unnecessary personnel.

6.1.1. For non-emergency personnel

Emergency procedures : Ventilate spillage area. Avoid contact with skin, eyes and clothing.

6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information

refer to section 8: "Exposure controls/personal protection".

Emergency procedures : Ventilate area. Stop leak if safe to do so.

6.2. Environmental precautions

Avoid release to the environment. Prevent liquid from entering sewers, watercourses, underground or low areas.

6.3. Methods and material for containment and cleaning up

For containment : Does not require any specific or particular technical measures.

Methods for cleaning up : Take up liquid spill into absorbent material.

Other information : Dispose of materials or solid residues at an authorized site.

6.4. Reference to other sections

For further information refer to section 13

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Keep container tightly closed. Ensure good ventilation of the work station. Wear personal

protective equipment. Avoid contact with eyes.

Hygiene measures : Do not eat, drink or smoke when using this product. Always wash hands after handling the

product

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.

7.3. Specific end use(s)

See Section 1.

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401



SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Glycerin (99.5%), chem. pure (56-81-5)		
United Kingdom - Occupational Exposure Limits		
Local name	Glycerol	
WEL TWA (OEL TWA) [1]	10 mg/m³ mist	
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE	

Monitoring methods		
Monitoring methods	A specific exposure sampling method is not available	
Biological monitoring methods	A specific exposure sampling method is not available	

8.2. Exposure controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

Personal protective equipment:

Data of item 8 do partly not refer to the use and the regular employing of the product (in this sense consult package leaflet and expert information), but to liberation of major amounts in case of accidents and irregularities.

Materials for protective clothing:

No specific measures are necessary

Hand protection:

None under normal use

Eye protection:

No specific measures are necessary

Skin and body protection:

No specific measures are necessary

Respiratory protection:

No specific measures are necessary

Environmental exposure controls:

Avoid release to the environment.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Viscous.

Colour : colourless to yellowish.

Odour : Odourless.
Odour threshold : No data available

pH : 4-9

Relative evaporation rate (butylacetate=1) : No data available

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401

Melting point : No data available Freezing point : No data available : No data available Boiling point Flash point : No data available Auto-ignition temperature : No data available Decomposition temperature : No data available Flammability (solid, gas) : Non flammable. Vapour pressure : No data available Relative vapour density at 20°C : No data available

Relative density : No data available Density : $\approx 1 \text{ g/cm}^3$ Solubility : Water: Miscible Log Pow : No data available Viscosity, kinematic : No data available

Viscosity, dynamic : 10000 – 25000 mPa·s 25°C Explosive properties : Product is not explosive.

Oxidising properties : Not oxidising.

Explosive limits : No data available

9.2. Other information

VOC content : 0 %

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

No materials to be especially mentioned.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral) : Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal) : Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation) : Not classified (Based on available data, the classification criteria are not met)

polyhexamethylene biguanide hydrochloride; PHMB (27083-27-8)	
LD50 oral rat	501 mg/kg
LD50 dermal rat	> 5000 mg/kg
LC50 Inhalation - Rat (Dust/Mist)	0.37 mg/l/4h

B BRAUN

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401

Skin corrosion/irritation

: Not classified (Based on available data, the classification criteria are not met)

B BRAUN

pH: 4 - 9

Serious eye damage/irritation : Not classified (Based on available data, the classification criteria are not met)

pH: 4 - 9

Respiratory or skin sensitisation : Not classified (Based on available data, the classification criteria are not met)

: Not classified (Based on available data, the classification criteria are not met) Germ cell mutagenicity

: Not classified (Based on available data, the classification criteria are not met) Carcinogenicity

Reproductive toxicity : Not classified (Based on available data, the classification criteria are not met)

STOT-single exposure : Not classified (Based on available data, the classification criteria are not met)

STOT-repeated exposure : Not classified (Based on available data, the classification criteria are not met)

Aspiration hazard : Not classified (Based on available data, the classification criteria are not met)

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : Harmful to aquatic life with long lasting effects.

Hazardous to the aquatic environment, short-term

(acute)

Hazardous to the aquatic environment, long-term

(chronic)

: Based on available data, the classification criteria are not met

: Harmful to aquatic life with long lasting effects.

polyhexamethylene biguanide hydrochloride; PHMB (27083-27-8)		
LC50 fish 1	0.026 mg/l Oncorhynchus mykiss, 96 h	
EC50 Daphnia 1	0.09 mg/l Daphnia magna, 48 h, (OECD 202 method)	
ErC50 algae	0.0191 mg/l Pseudokirchneriella subcapitata, 72h, (OECD 201 method)	
NOEC chronic crustacea	0.0084 mg/l Daphnia magna, 21 d, (OECD 211 method)	

12.2. Persistence and degradability

Prontosan Wound Gel X	
Persistence and degradability	The product has not been tested.

12.3. Bioaccumulative potential

Prontosan Wound Gel X	
Bioaccumulative potential	The product has not been tested.

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

Prontosan Wound Gel X

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

12.6. Other adverse effects

Other adverse effects : Slightly hazardous to water.

Additional information : Do not flush into surface water or sewer system

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401



SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Dispose of contents/container in accordance with licensed collector's sorting instructions.

Product/Packaging disposal recommendations

: Empty containers should be taken for recycling, recovery or waste in accordance with local regulation. Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse. Packaging that cannot be cleaned should be

disposed of like the product.

European List of Waste (LoW) code

07 06 99 - wastes not otherwise specified

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA / ADN / RID

ADR	IMDG	IATA	ADN	RID
14.1. UN number				
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
14.2. UN proper shippin	g name			
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
14.3. Transport hazard o	14.3. Transport hazard class(es)			
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
14.4. Packing group	14.4. Packing group			
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
14.5. Environmental hazards				
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
No supplementary information available				

14.6. Special precautions for user

Overland transport

Not regulated

Transport by sea

Not regulated

Air transport

Not regulated

Inland waterway transport

Not regulated

Rail transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

Contains no substance(s) listed on the REACH Candidate List

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

Safety Data Sheet

according to UK REACH Regulation SDS No: 00056-0401



Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

VOC content : 0%

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes:

All chapters have been modified since the previous version.

Abbreviations and acronyms:	
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
BLV	Biological limit value
BOD	Biochemical oxygen demand (BOD)
COD	Chemical oxygen demand (COD)
DMEL	Derived Minimal Effect level
DNEL	Derived-No Effect Level
EC-No.	European Community number
EC50	Median effective concentration
EN	European Standard
IARC	International Agency for Research on Cancer
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
NOAEL	No-Observed Adverse Effect Level
NOEC	No-Observed Effect Concentration
OECD	Organisation for Economic Co-operation and Development
OEL	Occupational Exposure Limit
PBT	Persistent Bioaccumulative Toxic
PNEC	Predicted No-Effect Concentration
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet

Safety Data Sheet

according to UK REACH Regulation

SDS No: 00056-0401



STP	Sewage treatment plant
ThOD	Theoretical oxygen demand (ThOD)
TLM	Median Tolerance Limit
VOC	Volatile Organic Compounds
CAS-No.	Chemical Abstract Service number
N.O.S.	Not Otherwise Specified
vPvB	Very Persistent and Very Bioaccumulative
ED	Endocrine disrupting properties
DOT	Department of Transport
TDG	Transportation of Dangerous Goods
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
GHS	Globally Harmonized System of Classification, Labelling and Packaging of Chemicals
IBC-Code	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
MARPOL 73/78	MARPOL 73/78: International Convention for the Prevention of Pollution From Ships
ADG	Transport of Australian Dangerous Goods

Other information

: Data of sections 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities. The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge. The delivery specifications are contained in the corresponding product sheet. This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

Full text of H- and EUH-statements:		
Acute Tox. 2 (Inhalation)	Acute toxicity (inhal.), Category 2	
Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4	
Aquatic Acute 1	Hazardous to the aquatic environment – Acute Hazard, Category 1	
Aquatic Chronic 1	Hazardous to the aquatic environment – Chronic Hazard, Category 1	
Aquatic Chronic 3	Hazardous to the aquatic environment – Chronic Hazard, Category 3	
Carc. 2	Carcinogenicity, Category 2	
EUH208	Contains polyhexamethylene biguanide hydrochloride; PHMB. May produce an allergic reaction.	
Eye Dam. 1	Serious eye damage/eye irritation, Category 1	
H302	Harmful if swallowed.	
H317	May cause an allergic skin reaction.	
H318	Causes serious eye damage.	
H330	Fatal if inhaled.	
H351	Suspected of causing cancer.	
H372	Causes damage to organs through prolonged or repeated exposure.	
H400	Very toxic to aquatic life.	
H410	Very toxic to aquatic life with long lasting effects.	
H412	Harmful to aquatic life with long lasting effects.	

Safety Data Sheet

according to UK REACH Regulation SDS No: 00056-0401



Skin Sens. 1B	Skin sensitisation, category 1B
STOT RE 1	Specific target organ toxicity – Repeated exposure, Category 1

Classification and procedure used to derive the classification for mixtures according to GB CLP Regulation		
Aquatic Chronic 3	H412	

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should therefore not be construed as guaranteeing any specific property of the product.