

PRODUCT DATA SHEET

Flexima® 3S Drainable Bags

Two-piece stoma bags for people with ileostomy or colostomy having pasty or liquid stools



Flexima® 3S is a 2-piece range with a guided mechanical coupling system for an easy positioning of the stoma bag onto the wafer. Its ergonomically shaped drainable bags have an integrated Roll'up soft clamp. They are available in two different sizes, midi and maxi.

Intuitive and secure closing system:

- Outlet with integrated self-gripping soft clamp
- 7 cm wide outlet opening
- The edges are softened to increase the wearing comfort once the Roll'up outlet is closed
- Once folded and closed, the drainable bag has a shape as compact as of a colostomy bag.

High performance protected filter (HP):

- Improved deodorization capacity through its high-tech breathable membrane
- Leak-proof and permeable to gases

Colors and capacity (volume as per ISO 8670-2):

- Midi:
 - Ø45 mm: capacity of +/- 650 ml, beige and transparent
 - Ø55 mm: capacity of +/- 635 ml, beige and transparent
 - Ø65 mm: capacity of +/- 615 ml, beige and transparent

- Maxi:
 - Ø45 mm: capacity of +/- 805 ml, beige and transparent
 - Ø55 mm: capacity of +/- 785 ml, beige and transparent
 - Ø65 mm: capacity of +/- 755 ml, beige and transparent
 - Ø80 mm: capacity of +/- 715 ml, beige and transparent

Material:

- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)

Others:

- To be used with the Flexima® 3S wafer corresponding to the chosen diameter of the bag
- Box of 30 bags

Advantages

Flexima® 3S stand for Simple, Soft and Secure:

- **Simple** coupling system that intuitively guides you through the attachment process so that only a minimum of patient education is required. Depending on which of the three slots you choose a multiple bag orientation is possible to adapt the bag positioning to your body contours.
- **Soft** coupling system with a floating flange that requires a minimum of pressure on the abdomen to prevent discomfort during application. The coupling system remains flexible when the wafer and bag are coupled together.
- **Secure** coupling system for three reasons: Enhanced leak-proof thanks to the close fitting of the bag on the wafer, increased security through the click when wafer and bag are well locked together and finally even more fixation with the possibility to use a belt attached to one of the three hooks on the bag.

Article Description	Article Number
FLEXIMA 3S ILEO MIDI BEIGE 45MM	932045N
FLEXIMA 3S ILEO MIDI BEIGE 55MM	932055N
FLEXIMA 3S ILEO MIDI BEIGE 65MM	932065N
FLEXIMA 3S ILEO MIDI TRP 45MM	932145N
FLEXIMA 3S ILEO MIDI TRP 55MM	932155N
FLEXIMA 3S ILEO MIDI TRP 65MM	932165N
FLEXIMA 3S ILEO MAXI BEIGE 45MM	932445N
FLEXIMA 3S ILEO MAXI BEIGE 55MM	932455N
FLEXIMA 3S ILEO MAXI BEIGE 65MM	932465N

Article Description	Article Number
FLEXIMA 3S ILEO MAXI BEIGE 80MM	932480N
FLEXIMA 3S ILEO MAXI TRP 45MM	932545N
FLEXIMA 3S ILEO MAXI TRP 55MM	932555N
FLEXIMA 3S ILEO MAXI TRP 65MM	932565N
FLEXIMA 3S ILEO MAXI TRP 80MM	932580N