Venous Access Port System




Vascular Systems




Celsite[®] Safety

Access ports for mid to long-term venous applications



PEEK and titanium combination as special safety features



Long Term Central Venous Access without Compromise on Safety

Large puncture area with high density silicone septum

For simple puncturing and reliable sealing to allow good port life

Extra large suture holes

For easy fixation of the access port with sutures

Lasermarked PUR and silicone catheters

- Clear readability of the catheter length
- No ink is added

Radiopaque connection ring

Anti kink protection and additional fixation of the catheter



High pressure resistant

- Complete range of Celsite® Safety is resistant to high pressure injections up to 325 psi
- Enables power injections of contrast media
- No need for additional venous access



Radiopaque CT - marking

Clear identification of high pressure resistance under the x-ray





Celsite[®] Safety

Long Term Central Venous Access without Compromise on Safety



Surecan® Safety II - Port Needle

- Intuitive safety mechanism to reduce the risk of needle stick injuries
- High pressure resistance up to 325 psi






Safecan™ Safety - Puncture Needle

- Safety mechanism to reduce the risk of needle stick injuries
- Echogenic puncture needle for needle tip location via ultrasound



Accessories



	Implantation technique	Seldinger	Surgical cut-down	Separate Order Numbers		
	Kit designation	1	2			
Α	Tunnelling Rod	✓	-	-		
В	Tear-away Introducer	L 180/140 mm	-	-		
С	Omnifix® Luer Syringe	10 mL	-	4616103V		
D	Straight Surecan®	22 G x 30 mm	22G x 30 mm	4439848		
E	Vein lifter	✓	~	-		
F	Rinsing Hub	V	~			
G	Safecan™ Safety - Puncture Needle	18 G x 70 mm	-	613104		
Н	Surecan® Safety II	20 G x 20 mm	-	4447006		
1	J Guide Wire with Dispenser	0.035" x 50 cm	-	-		

Celsite® Safety offers a wide range of Silicone and PUR catheters as well as two different port sizes, Standard and Small.

Catheter	OD (F/mm)	ID (mm)	Length (mm)	Flow rate* (ml/min)		Recommended maximum flow rates (mL/s) PSI Contrast media at 37°C (325 psi = 22.4 bar)**			Implantation technique	Туре	Reference	Accesso- ries see page 3
	(F/IIIII)					Viscosity up to 11.4 mPa.s (cP)						
				19 G	22 G	22 G	20 G	19 G				
Standard												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Surgical cut-down	T601F	4437556	2
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Surgical cut-down	T601L	4437573	2
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Surgical cut-down	T601P	4437565	2
PUR	8.5 / 2.8	1.6	500	48	12	2	5	5	Surgical cut-down	T601H	4437581	2
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SST601F	4437603	1
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SST601L	4437612	1
Silicone	10 / 3.2	1.6	500	48	12	2	5	5	Seldinger	SST601G	4437620	1
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Seldinger	SST601P	4437607	1
PUR	8.5 / 2.8	1.6	500	48	12	2	5	5	Seldinger	SST601H	4437617	1
Small												
Silicone	6.5/2.2	1.1	500	26	10	2	5	5	Surgical cut-down	T605F	4437758	2
Silicone	6.5/2.2	1.1	500	26	10	2	5	5	Seldinger	SST605F	4437803	1
Silicone	8.5/2.8	1.2	500	34	11	2	5	5	Seldinger	SST605L	4437817	1
PUR	5/1.6	1.1	500	26	10	2	5	5	Seldinger	SST605C	4437800	1
PUR	6.5/2.1	1.4	500	37	12	2	5	5	Seldinger	SST605P	4437809	1
PUR	8.5/2.8	1.6	500	48	12	2	5	5	Seldinger	SST605H	4437813	1


^{*} Gravity infusion of saline (0.9%) through a 22G respectively 19G needle from a height difference of 1 m and a catheter length of 40 cm. According to ISO 10555-1

All Celsite® Safety Ports are PVC, LATEX and DEHP free









^{**} Flow rates determined according to ISO 10555-6 with a catheter of 20 cm and Surecan® Safety II 20 mm needle


Customised Add-on kits for Access Ports



- B. Braun offers individual add-on kits which contain an organized array of disposable components for your specific port implantation procedure.
- B. Braun customised add-on kits save procedural time since they facilitate a quick and easy setup and cleanup.
- B. Braun customised add-on kits are cost effective by simplifying your purchasing and accounting procedures and by consolidating components in one affordable kit.

B. Braun customised add-on kits mean:

- Custom made kit solutions
- All components in one kit
- Less preparation times
- More time for the patient
- Compact storage size
- Reduced waste volume
- Known and same quality
- One supplier contact



Distributor

B. Braun Melsungen AG | Vascular Systems | Sieversufer 8 | 12359 Berlin | Germany Phone +49 30 568207-300 | Fax +49 30 568207-130 | www.bbraun.com

Manufacturer acc. to MDD 93/42

B. Braun Médical \mid 204, avenue du Maréchal Juin \mid 92107 Boulogne Cedex \mid France Tél +33 1 41 10 53 00 \mid Fax +33 1 41 10 52 05 \mid www.bbraun.fr

Aesculap AG \mid Am Aesculap-Platz \mid 78532 Tuttlingen \mid Germany Phone +49 74 61 95-0 \mid Fax +49 74 61 95-26 00 \mid www.aesculap.com

Aesculap - a B. Braun company

The main product trademark 'Aesculap'is a registered trademark of Aesculap AG.
The product trademarks "Celsite", "Surecan" and "Cytocan" are registered trademarks of B. Braun Melsungen AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure no. 6050340

1215/1.5/1