



ACTREEN® MINI RANGE

SMALL IS BEAUTIFUL... Female self intermittent catheters

The new Mini range

A woman's life should be rich and exciting, full of personal & professional hopes and aspirations.

Every woman has her own plan for fulfilling these dreams, but it's often the small things that make a big difference to their daily lives.

The "Mini" range from B. Braun has been designed to make life easier, every day, everywhere, offering the freedom and peace of mind you are looking for with simple and convenient to use products.

So, whatever type of woman you are we will help you find the most suitable catheter!

Whether you are looking for discretion, style, easeof-use or safety, or whether you are able-bodied or disabled, suffering from urinary disorders with a specific need for urinary catheters, we have the answer...

Actreen® Mini range is your perfect indoor and outdoor companion, always ready when you are...





Are you Mini Set...



New innovative set

All in one product: Catheter + bag & non return valve + discreet pouch

Lightest Mini set weight*





SIMPLIFIED PROCEDURE

- 1 Unfold the set
- 2 Open the catheter
- 3 Start catheterization



Discreet pouch

MINI "PHILOSOPHY"

- Straight forward use, no need for education
- Ready to use, no need for water, no mess, no splash!
- Low profile products, discreet into a pocket, in a handbag
- Light weight, easy to carry
- Optimal waste management in weight & volume*
- A discreet pouch into each box!

...or Mini Cath ?



Mini Cath has a new design, more feminine look, stylish + discreet pouch...









- 1 Open the catheter
- 2 Start catheterization



Discreet pouch






The Mini Range

Common benefits...

The Mini catheters benefit from the same technology and design:

- Non traumatic eyes highly positioned for effective drainage
- Safe catheter material, PVC free, DEHP free
- Ready to use catheter, pre-lubricated with Hydrophilic lubricant
- Guarantees the same catheter hydration for over 60 min from opening*!
- No need for water, no liquid management, no waiting time!



... in 2 unique concepts

Actreen® Mini Cath

Actreen® MiniSet





	Actreen® Mini Cath	REFERENCE	UNITS/BOX
And I was	■ 10 (3,3 mm)	228010*	30 +1 discreet pouch
	☐ 12 (4,0 mm)	228012*	30 +1 discreet pouch
	■ 14 (4,7 mm)	228014*	30 +1 discreet pouch
	Actreen® Mini Set	REFERENCE	UNITS/BOX
	Actreen® Mini Set ■ 10 (3,3 mm)	239010*	UNITS/BOX 30 +1 discreet pouch
Tan Te			•

Manufactured and distributed by:

B. Braun Medical | 204 avenue du Maréchal Juin | 92100 Boulogne-Billancourt | FRANCE

Tél. 01 41 10 53 00 | Fax 01 41 10 53 99 | www.bbraun.com

B. Braun Medical | SAS with a share capital of 31 000 000 € | RCS Nanterre 562050856

ZG00730 - 04/2016